

CLARITY – LIST OF CURRENTLY CONTROLLED INSTRUMENTS (version 8.3)

GCs (p/n A23)

	Producer	Name	Interface	Status
	Agilent	4890D, 5890 II, 5890A, 6890	RS232	Ready
	Agilent	6820, 6850, 6850 II, 6890 Plus, 6890N	RS232 or LAN	Ready
New	Agilent	6850, 6890 Plus, 6890A, 6890N, 7820, 7890A, 7890B, 8860, 8890, Intuvo 9000	LAN - ICF	Testing
	Agilent	7890A, 7890B	LAN	Ready
	Ample Technology Center	ATC-6900 GC	LAN	Ready
	Apix ^{/1,2}	ChromPix2, TwinPix	LAN	Testing
	Dani	GC1000	RS232 + A/D converter	Ready
	Dani ^{/1,3}	Master GC	RS232 or LAN	Ready
	Ellutia ^{/1,2}	200 Series	RS232 + A/D converter	Ready
	Ellutia ^{/1,3}	300 Series	RS232 + A/D converter	Ready
New	Full ^{/1}	GC 9720 Plus, GC 9790 Plus	LAN	Ready
	GOW-MAC	Series 816	LAN	Ready
	NETEL ^{/1}	Analyte 2900A, Chrom Lite 3000A	A/D converter + RS232	Testing
	Shimadzu	GC14C, GC17A	RS232, OPT-USB + A/D converter	Ready
New	Shimadzu	GC2010, GC2010 Plus, GC2010 Pro, GC2014, GC2014C, GC2014C APC/AFC	RS232	Ready
	Snir ^{/1,2}	Sion 4210	LAN	Testing
New	YoungIn Chromass ^{/1}	ChroZen GC	LAN	Testing
	YoungIn Chromass ^{/1}	YL6100	RS232	Ready
	YoungIn Chromass ^{/1}	YL6500	LAN	Ready

HPLC Systems (p/n A24)

	Producer	Name	Interface	Status
	Agilent	1100 LC, 1200 LC, 1260 Infinity, 1290 Infinity II (see list of particular components)	LAN - ICF	Ready
	Agilent	1100 system	GPIB, LAN	Ready
	Agilent	1120 LC (G4286A, G4287A, G4288A, G4289A, G4290A), 1220 LC (G4286B, G4286C, G4287B, G4287C, G4288B, G4289B, G4290B, G4291B, G4291C, G4292B, G4292C, G4293B), 1220 VL (G4288C, G4289C, G4290C, G4293C, G4294B)	LAN - ICF	Ready
	Agilent	1200 system (selected modules only)	LAN	Ready
	Agilent	1260 Infinity II	LAN - ICF	Testing
	Alltesta ^{/1}	Analyzer	RS232	Testing
	Analytical Technologies Limited ^{/4}	iUHPLC Series 3000Plus	RS232	Testing
	Antec	ALEXYS	RS232	Ready
	Antec ^{/5}	ROXY	RS232	Ready

	CoruiTech ^{/1}	Rainbow, RainbowC	LAN	Testing
	CVC Technologies ^{/4}	Atima	LAN	Testing
	CVC Technologies ^{/4,6}	Laevitas	USB	Ready
	Dalian Elite ^{/1,7}	EClassical 3100, iChrom 5100, SOPTOP 1620A Plus	LAN	Ready
New	Esensing ^{/4}	MassChrom 1299	RS232 + LAN	Testing
	GL Science ^{/1,6}	LC800	USB	Ready
	Hanon Instruments ^{/4}	LC7000, LC7300	LAN/RS232	Ready
	Hitachi ^{/8,9}	Chromaster	USB	Ready
	Hitachi ^{/6,9}	LaChrom Elite	USB	Ready
	Hitachi ^{/6,8,9}	Primaide	USB	Ready
	Ingos ^{/1}	LC 5000	LAN	Testing
	Interlab ^{/2,9}	MAESTRO	USB	Ready
	Knauer ^{/1}	Azura, BlueShadow, Smartline, WellChrom	RS232	Ready
	PerkinElmer ^{/1}	Flexar	RS232 + LAN	Testing
	PerkinElmer	NexSAR	RS232	Testing
	RIGOL ^{/1}	L-3000	RS232	Testing
	Shimadzu ^{/6}	LC-10, LC-20 Prominence	RS232 or LAN + corresp. controller (see HPLC Other)	Ready
New	Spark	LC1299, LC1299 Plus	RS232, USB or LAN	Testing
	SSI	B2300, VersaGrad, VersaGrad UHP	RS232 + LC UNI	Ready
	Sykam ^{/1}	S500, S600	RS232	Ready
	YoungIn Chromass ^{/1}	ChroZen	LAN	Ready
	YoungIn Chromass ^{/1}	YL9100, YL9300	LAN/RS232	Ready
	YoungIn Chromass ^{/4}	YL9800	USB	Ready

HPLC Pumps (p/n A24)

	Producer	Name	Interface	Status
	Agilent	1100 Binary, 1100 Isocratic, 1100 Quaternary	GPIB, LAN	Ready
	Agilent	1200 Binary, 1200 Isocratic, 1200 Quaternary	LAN	Ready
	Agilent	Binary (G1312A, G1312B, G1312C), Isocratic (G1310A, G1310B), Quaternary (G1311A, G1311B, G1311C, G4204A) (all including selected pump valve clusters)	LAN - ICF	Ready
New	Agilent	Binary (G1376A, G2226A, G4220A, G4220B, G4302A, G4782A, G7112B, G7120A, G7161A, G7161B), Isocratic (G1361A, G7110B), Quaternary (G5611A, G5654A, G7104A, G7104C, G7111A, G7111B), u G1776A (all including selected pump valve clusters)	LAN - ICF	Testing
	Alltech/Grace	301, 426, 626, 726	RS232 - LC UNI	Ready
	Analytical Technologies Limited	ALC2010, ALC2100, COLUMN-WASH	RS232 - LC UNI	Ready
	Analytical Technologies Limited ^{/4}	P-3210U, PB-3220U, PQ-3245U	RS232	Testing
	Antec	LC100	RS232	Ready
	Antec ^{/4}	LC110	RS232	Ready
New	Antec ^{/4}	LC110S	RS232 or LAN	Ready
	ASI	500, 501, 520	RS232 - UNI Ruby	Ready

	Beckman/Altex	110A, 110B, 114	Zebrick	Ready
	Bio-Rad	1350	Zebrick	Ready
	Bischoff	2250	RS232	Ready
	Chrom Tech	ISO-100, ISO-1000, ISO-1500, ISO-2000	RS232	Ready
	Chromophor	Chromopeak 2010	RS232 - UNI Ruby	Ready
	Chromsystems	CLC 300	Zebrick	Ready
	Chromsystems	CLC 330	RS232	Ready
New	CQS ^{/4}	Climax 150 Pump, Climax 155 Pump, Climax S 1130	RS232	Ready
	Ecom ^{/1}	Alpha 05, Alpha 10, Alpha 100, Alpha 250, Alpha 50, Beta 10, Beta 50, IOTA-250, IOTA-300, IOTA-400, IOTA-S-10, IOTA-S-50, IOTA-S-100, IOTA-S-200, IOTA-S-300, KAPPA10PP, LCP 4050, LCP 4100	RS232	Ready
	Ecom	Alpha 10, Alpha 100, Alpha 50, Beta 10, Beta 50, KAPPA10, LCP 4100	RS232 - LC UNI	Ready
	Ecom ^{/1}	ECP2002, ECP2010, ECP2050, ECP2100, ECP2300	RS232 or LAN	Ready
	Ecom	LCP 4000, LCP 4020	Zebrick	Ready
New	Esensing	SPH1299	RS232, USB or LAN	Testing
	Flom	301M, AI-12 series	RS232 - UNI Ruby	Ready
	Flom	UI-22 (UI series)	RS232	Testing
	Flux - Thermo ^{/1}	Allegro, Rheos 2000, Rheos 2200	RS232	Ready
New	Futecs ^{/1}	P-6000	RS232 - UNI Ruby	Testing
	Gilson	302, 303, 305, 306, 307	GSIOC	Ready
	Gynkotek	M300CS, M480	Zebrick	Ready
	Hanon Instruments ^{/4}	LC7010, LC7011, LC7012	LAN	Ready
	Hitachi ^{/8,9}	Chromaster 5110, 5160	USB	Ready
	Hitachi ^{/9}	L-2100, L-2130	USB	Ready
	Hitachi	L-6000	via L-6200	Ready
	Hitachi	L-6200, L-6210	RS232	Ready
	Hitachi ^{/6,8,9}	Primaide 1110	USB	Ready
	JASCO	PU-880	Zebrick	Ready
	Knauer ^{/1}	10P/20P, 40P, 80P, P2.1S/P4.1S, PLATINblue P-1, S100, Smartline 1050	RS232 or LAN	Ready
	Knauer ^{/1}	K-1000, K-1001, K-120, K-1800, Prep Pump 1800, K-501, Smartline S1000	RS232	Ready
	Knauer	K-500, K-501, M-64	Zebrick	Ready
	Knauer ^{/1}	P2.1L, P6.1L	LAN	Ready
	Knauer	Smartline S1000	RS232 - LC UNI	Ready
	Kontron	320, 322, 325	RS232 - LC UNI	Testing
	Kontron	420, 422, 525	RS232 - LC UNI	Ready
	LKB - Amersham	LKB 3500, P-50, P-500	Zebrick	Ready
	PerkinElmer ^{/1}	Flexar Binary LC Pump, Flexar Quaternary LC Pump	RS232	Testing
	PerkinElmer	NexSAR pump	RS232	Ready
	Postnova Analytics	PN1130	Zebrick	Ready
	Recipe	GP3000, IP3000	RS232	Ready
	RIGOL ^{/1}	L-3200	RS232	Testing
	Schambeck	9404, SDS 9414 I	Zebrick	Ready
	Schambeck ^{/4}	S 9425	RS232	Ready
	Schambeck	SFD 2100	RS232	Ready

Shimadzu	LC-10AD, LC-10ADvp, LC-10Ai, LC-10AS, LC-10AT, LC-10ATvp, LC-20AB, LC-20AD, LC-20ADXR, LC-20AS, LC-20AT, LC-6AD, LC-8A	Shimadzu controller (see HPLC Other)	Ready
Shimadzu	LC-10ADvp, LC-10ATvp	OPT-USB	Ready
Shimadzu	LC-20AB, LC-20AD, LC-20ADXR, LC-20AS, LC-20AT (in LC-10 emulation mode only)	OPT-USB	Ready
Shimadzu	LC-6A	Zebrick	Ready
Spark	SPH1240	RS232, USB or LAN	Ready
Spark	SPH1299	RS232, USB or LAN	Testing
SSI	150, 1500, HF 300, Prep 100, Prep 300, Q-Grad, Series I, Series II, Series III, UHP, VersaGrad UHP	RS232 - UNI Ruby	Ready
SSI	222D	Zebrick	Ready
SSI	LD-Class, LS-Class, LU-Class, M1-Class, MX-Class, PR-Class	RS232	Ready
SSI	Mighty-Mini	RS232 - UNI Ruby	Testing
Sykam	S1021, S1521, S1122	Zebrick	Ready
Sykam ^{/1}	S1125, S1130, S1132	RS232	Ready
Sykam	S2100	RS232	Ready
Techlab	Minipump1	RS232 - UNI Ruby	Ready
Teledyne ISCO	2300	Zebrick	Ready
TSP	Constametric III	Zebrick	Ready
Unicam	PU 4015	Zebrick	Ready
Varian	212-LC	RS232 - UNI Ruby	Ready
Waters	6000A, M 45, 501, 510, 515	Zebrick	Ready
Watrex	DeltaChrom P200, DeltaChrom SDS20, DeltaChrom SDS30	RS232 - LC UNI	Ready
Watrex	Mighty-Mini	RS232 - UNI Ruby	Testing
New Watrex ^{/1}	Streamline P1m, Streamline P1s	USB	Testing
YoungIn Chromass	ChroZen UHPLC Pump	LAN	Testing
YoungIn Chromass ^{/1}	YL9110, YL9111, YL9112	LAN	Ready
	Any pump with flow controlled by voltage or frequency	CB20 or Zebrick	Ready
	Any pump with on-line flow controlled by text (ASCII) comm. protocol	RS232 - LC UNI	Ready

HPLC Detectors (p/n A24)

Producer	Name	Interface	Status
Agilent	1100 MWD, 1100 RID, 1100 VWD	GPIB, LAN	Ready
Agilent ^{/10}	1200 MWD	LAN	Ready
Agilent	1200 VWD	LAN	Ready
Agilent	1200: ELSD	RS232	Ready
Agilent	ELSD (G4218A, G4260A, G4260B, G4261A, G4261B, G7102A)	LAN - ICF	Testing
Agilent	FLD (G1321A, G1321B)	LAN - ICF	Ready
Agilent	FLD (G1321C, G7121A, G7121B)	LAN - ICF	Testing
Agilent	MWD (G1365A, G1365B, G1365C, G1365D), RID (G1362A), VWD (G1314A, G1314B, G1314C, G1314D)	LAN - ICF	Ready
Agilent	MWD (G7165A), RID (G7162A, G7162B), VWD (G1314E, G1314F, G7114A, G7114B)	LAN - ICF	Testing
Alltech/Grace	3300 ELSD	LAN	Ready

	Analytical Technologies Limited ^{/4}	ASPD2600	RS232	Ready
	Analytical Technologies Limited ^{/4}	UVD-3500U	RS232	Testing
	ANALYTIK Jena ^{/1}	PlasmaQuant MS	LAN	Testing
New	ANALYTIK Jena ^{/4}	PQ LC IC System 150	RS232	Testing
	Antec	DECADE Elite, DECADE Lite	USB or LAN	Ready
	Antec	DECADE II, DECADE II MD, DECADE II SDC	RS232	Ready
	Antec ^{/5}	ROXY Potentiostat	RS232	Ready
New	Chromophor ^{/4}	C 150 Ion Chromatography System	RS232	Testing
	CoruiTech ^{/1}	Rainbow detector	USB	Testing
New	CQS ^{/4}	Climax 150 IC, Climax 155 IC	RS232	Testing
New	CQS ^{/4}	Climax 150 UV, Climax S 3250 UV/Vis detector	RS232	Ready
	Ecom ^{/1}	Apatite, Flash 06 DAD, Flash 06 Single, Flash 10 DAD, Malachite, Opal, Sapphire	RS232	Ready
	Ecom ^{/1}	BABY18DAD, TOY14DAD, TOY18DAD	RS232 or LAN	Testing
	Ecom ^{/1}	ECD2600, ECD2800	RS232 or LAN	Ready
	Ecom	RI2012	RS232 - UNI Ruby	Ready
	Ecom	Topaz	RS232	Ready
	Ecom	TOYFIX	RS232	Testing
	Elysia-raytest ^{/1,11}	Gabi Nova	USB or LAN	Testing
	Elysia-raytest ^{/1,11}	Gabi Star Radiodetector	RS232	Ready
	Elysia-raytest ^{/1,11}	POMO Radiodetector	USB - UNI Ruby	Ready
	Elysia-raytest ^{/1,11}	Ramona Star	RS232	Testing
	ERC	RefractoMax 520 Series	RS232 - UNI Ruby	Ready
	ESA	Coulochem III	USB	Ready
New	Esensing ^{/4}	ELSD	RS232	Testing
New	Esensing ^{/4}	MCH UV	LAN	Ready
New	Futecs ^{/1}	CD-6000, RI-6000	RS232 - UNI Ruby	Testing
	Gilson	PrepELS™ II	RS232	Ready
	Hanon Instruments ^{/4}	LC7020	LAN	Ready
	Hanon Instruments ^{/4}	LC7070, LC7081	RS232	Ready
	Hi Corporation	RI2012	RS232 - UNI Ruby	Ready
	Hitachi ^{/8,9}	Chromaster 5410, 5420, 5440, 5450	USB	Ready
	Hitachi ^{/9}	L-2400	USB	Ready
	Hitachi ^{/6,8,9}	Primaide 1410	USB	Ready
	Interlab ^{/2}	MAESTRO ELSD	RS232	Testing
	Interlab	Omega FLD, Omega UV-VIS	USB	Testing
	Japan Analytical Industry Co., LTD. ^{/1}	RI-700II	RS232	Ready
New	JASCO	FP-4020, FP-4025	RS232	Testing
	Knauer ^{/1}	10D, 40D, 50D, CM2.1S, PLATINblue MW-1, S2520 VWD, S2550 VWD, UVD2.1L, UVD2.1S	RS232 or LAN	Ready

	Knauer ^{/1}	Conductivity Monitor 2900, K-2600, S200 UV, S2300 RID, S2400 RID, Smartline S2500 VWD, WellChrom RI (K-2300, K-2301, K-2400, K-2401), WellChrom UV (K-2000, K-2001, K-2500, K-2501)	RS232	Ready
	Knauer ^{/1}	MWD2.1L	LAN	Ready
	Knauer	RID2.1L	LAN	Testing
	Knauer ^{/1,12}	S2510	RS232	Ready
	Kontron	432	A/D converter + RS232 - UNI Ruby	Ready
	PerkinElmer ^{/1}	Flexar FX-UV/Vis, Flexar UV/Vis	RS232	Testing
	PerkinElmer ^{/1}	NexION 1000, NexION 2000	LAN	Testing
	Recipe	EC6000	RS232 - UNI Ruby	Testing
	Recipe ^{/4}	FW3000, UV3000	RS232	Ready
	RIGOL ^{/1}	L-3500	RS232	Testing
	RIGOL	L-3530	RS232	Ready
	Schambeck	RI 2000, RI 2012, S 2020, S 2030	RS232 - UNI Ruby	Ready
	Schambeck ^{/4}	S 4245, ZAM 3000, ZAM 4000	RS232	Ready
	Schambeck	S3210, S3240	RS232	Testing
	Sedere	Sedex 100LT, Sedex 90LT, Sedex FP, Sedex LC	RS232 or USB	Ready
	Sedere ^{/13}	Sedex 80LT	RS232	Ready
	Sedere	Sedex 85LT	RS232	Ready
	Semba Biosciences ^{/4}	Octave 4X, Octave 4XE	RS232	Ready
	Shimadzu	ELSD-LTII	RS232	Ready
	Shimadzu	RF-10A, RF-10Axl, RF-20A, RF-20Axs, RID-10A, SPD-10A, SPD-10Ai, SPD-10AV, SPD-10AVi, SPD-10Avp, SPD-10AVvp, SPD-20A, SPD-20AV	Shimadzu controller (see HPLC Other)	Ready
	Shimadzu ^{/4}	RF-10Axl, RF-20Axs	RS232	Ready
	Shimadzu ^{/4,14}	RF-20A	RS232	Ready
	Shodex	RI-101, RI-102, RI-104	RS232 - UNI Ruby	Ready
	Shodex	RI-501, RI-501EX, RI-502, RI-504	LAN - UNI Ruby	Testing
	Softa	1300, 1400, 200S, 300S, 400	RS232	Ready
	Sykam	RI2000, S3585, S3590	RS232 - UNI Ruby	Ready
New	Sykam ^{/1}	S 150	RS232	Testing
	Sykam	S3210, S3240	RS232	Testing
	Sykam ^{/1}	S3245, S3250	RS232	Ready
	Unimicro ^{/1}	UNIEX-7700 ELSD	RS232	Testing
New	VICI - Valco ^{/1}	TCD 3	RS232	Ready
New	Watrex ^{/4}	Streamline IC1	USB	Testing
New	Watrex ^{/4}	Streamline UV1	RS232	Ready
	YoungIn Chromass ^{/1}	ChroZen HPLC-A UVD, ChroZen UHPLC UVD, YL9120	LAN	Ready
	YoungIn Chromass ^{/1}	YL9170, YL9180	RS232	Ready
	YoungIn Chromass	YL9181	RS232	Ready

PDA Detectors (p/n A24 + p/n A29)

Producer	Name	Interface	Status
Agilent	1100: DAD, FLD	GPIB, LAN	Ready
Agilent ^{/10}	1200: DAD, FLD	LAN	Ready
Agilent	FLD (G1321A, G1321B)	LAN - ICF	Ready
Agilent	FLD (G1321C, G7121A, G7121B)	LAN - ICF	Testing
Agilent	G1315A, G1315B, G1315C, G1315D, G4212A, G4212B	LAN - ICF	Ready
Agilent	G7115A, G7117A, G7117B, G7117C	LAN - ICF	Testing
Analytical Technologies Limited ^{/4}	DAD-3520U	RS232	Testing
Dalian Elite ^{/1,7}	iChrom 5100 PDA	LAN	Testing
Duratec ^{/15}	DDT-3200 PDA	USB	Ready
Ecom ^{/1}	ECDA2000, ECDA2800	USB or LAN	Testing
Ecom ^{/1}	Flash 12 DAD, FLASH 14 DAD	RS232	Ready
Ecom	TOYDAD 400, TOYDAD 400 L, TOYDAD 600, TOYDAD 600 L, TOYDAD 800, TOYDAD 800 L	RS232	Testing
<i>New</i> Esensing ^{/4}	MCH PDA	LAN	Ready
Hanon Instruments ^{/4}	LC7060	LAN	Ready
Hitachi ^{/8,9}	Chromaster 5430	USB	Ready
Hitachi ^{/6,8,9}	Primaide 1430	USB	Ready
Knauer ^{/1}	DAD2.1L, DAD6.1L, K-2700, K-2800, PLATINblue PDA -1, S2800 PDA, S2850 PDA	LAN	Ready
Knauer ^{/1}	Smartline S2600 MWD	RS232	Ready
PerkinElmer ^{/1}	Flexar Plus PDA	LAN	Testing
Prince ^{/1}	Next DAD	RS232	Testing
RIGOL ^{/1}	L-3520	USB or LAN	Testing
Schambeck	S3210, S3240	RS232	Testing
Shimadzu	SPD-M10Avp, SPD-M20A	Shimadzu controller (see HPLC Other)	Ready
Sykam	S3210, S3240	RS232	Testing
YoungIn Chromass ^{/1}	ChroZen HPLC-A PDA, ChroZen UHPLC PDA, YL9160	LAN	Ready

LC-MS Detectors (p/n A38 + p/n A24)

Producer	Name	Interface	Status
Advion	Expression CMS	USB	Ready
Microsaic	4000 MiD, 4500 MiD	LAN	Ready
YoungIn Chromass	YL9900	USB	Ready

GC-MS Detectors (p/n A37 + p/n A23)

Producer	Name	Interface	Status
YoungIn Chromass	YL6900	LAN	Testing

GC-MS TOF Detectors (p/n A37 + p/n A23)

Producer	Name	Interface	Status
Dani ^{/1}	MasterTOF	USB	Ready
Markes International	BenchTOF	LAN	Ready

HPLC Other (p/n A24)

Producer	Name	Interface	Status
Bronkhorst ^{/4}	mini Cori-Flow	RS232	Testing
ECOM	Panda-30 HID	USB	Ready
Knauer ^{/1}	ASM2.1L	LAN	Ready
Knauer ^{/1}	I-Box, IF2, Manager 40M, Smartline Manager 5000	RS232	Ready
Knauer ^{/1}	IFU 2.1	USB	Ready
Shimadzu	CBM-20A, CBM-20Alite, SCL-10Avp	RS232 or LAN	Ready
SISw	SSI pumps USB interface	USB - UNI Ruby	Ready
SISw ^{/1}	Systec Degasser Watchdog, USB pH monitor	USB - UNI Ruby	Ready
Unimicro ^{/1}	HVPS	RS232	Testing

GC Accessories (p/n A23)

Producer	Name	Interface	Status
Faculty of Science, Charles University in Prague	GCxGC Modulator	RS232	Ready

Autosamplers (p/n A26)

Producer	Name	Interface	Status
Agilent	1100 ALS, 1100 TAS	GPIB, LAN	Ready
Agilent	1200 ALS, 1200 TAS	LAN	Ready
Agilent	6850 ALS	through GC6850	Ready
New Agilent	7650, 7673B, 7673C, 7683A, 7683B, 7693A, ALS 6850	LAN - ICF	Testing
Agilent	7673, 7673 II	RS232	Ready
Agilent	7683	through GC6890	Ready
Agilent	7683, 7693	through GC7890	Ready
Agilent	G1313A, G1329A, G1329B (all including thermostat (G1330A, G1330B))	LAN - ICF	Ready
Agilent	G1327A, G1367A, G1367B, G1367C, G1367D, G1367E, G1368A, G1377A, G1378A, G1389A, G2258A, G2260A, G4226A, G4227A, G4303A, G4767A, G5667A, G5668A, G7129A, G7129B, G7129C, G7157A, G7167A, G7167B (all including thermostat (G1330A, G1330B))	LAN - ICF	Testing
Alltech/Grace	570, 580	RS232	Ready
Ample Technology Center	ATC3000A, ATC3100A, ATC3200A	LAN	Ready

	Analytical Technologies Limited ^{/4}	AS-3320U	RS232	Testing
	Analytical Technologies Limited	AS2000	RS232	Ready
New	ANALYTIK Jena ^{/4}	PQ LC 5250, PQ LC 5300	RS232	Ready
	Antec	AS 100	RS232	Ready
	Antec ^{/4}	AS 110	RS232 or LAN	Ready
	Antec ^{/4}	AS 110 Micro	RS232	Ready
	Cetac	ASX-8000	USB	Ready
New	Chromsystems	CLC 200, CLC 220	RS232	Ready
New	CQS ^{/4}	Climax 150 AS, Climax S 5250, Climax S 5300	RS232	Ready
	CTC Analytics	A200S	RS232	Ready
	CTC Analytics	PAL (GC, LC, Combi), PAL-xt (LC-xt, HTC-xt, HTS-xt, HTX-xt, GC-xt, Combi-xt, Combi-xt extended)	RS232 or LAN	Ready
	CTC Analytics	PAL LSI, PAL RSI, PAL RTC, PAL3 System	LAN	Testing
	Dani	ALS1000	RS232	Ready
	Dani ^{/1,3}	Master AS, Master DHS	RS232	Ready
	Dani ^{/1,3}	Master SHS	RS232, USB or LAN	Ready
	Eicom ^{/2}	AS-700	RS232, USB or LAN	Ready
	ESA	542	RS232	Ready
New	Esensing	MCH Alias	RS232, USB or LAN	Ready
	Hanon Instruments ^{/4}	LC7050	RS232	Ready
	Hitachi ^{/8,9}	Chromaster 5210, Chromaster 5260, Chromaster 5280	USB	Ready
	Hitachi ^{/9}	L-2200	USB	Ready
	Hitachi ^{/6,8,9}	Primaide 1210	USB	Ready
	HTA ^{/1}	HT1500L	USB	Ready
	HTA	HT2000H, HT2000HT, HT2100H, HT2800T, HT3000A, HT3100A, HT3200A, HT4000L	LAN	Ready
	HTA	HT200H, HT250D, HT280T (without SPME), HT300A, HT310A, HT300LV, HT310LV, HT300L, HT310L	RS232	Ready
	HTA ^{/1}	HT800L	RS232	Ready
	Knauer	AS 3800, AS 3900, Optimas	RS232	Ready
	Knauer ^{/1}	AS 3950	RS232	Ready
	Knauer	AS 6.1L	RS232 or LAN	Ready
	Knauer ^{/1}	PLATINblue AS-1	RS232 or LAN	Ready
	Kontron	360, 460, 465, 560	RS232 - UNI Ruby	Ready
	PerkinElmer ^{/1}	Flexar AS	RS232, USB, LAN	Testing
	PerkinElmer	NexSAR Inert HPLC Autosampler	RS232, USB or LAN	Ready
	Postnova Analytics	PN5300	RS232	Ready
	Recipe	AS3000, AS5000	RS232	Ready
	RIGOL ^{/1}	L-3320	RS232, USB or LAN	Testing
	Schambeck ^{/4}	S 6250, S 6300	RS232	Ready
	Schambeck	S5200	RS232	Ready
	Shimadzu	AOC-20i, AOC-20s, AOC-20d	RS232	Ready
New	Shimadzu	AOC-5000 Plus	RS232 or LAN	Ready

	Shimadzu	SIL-10ADvp, SIL-10AF, SIL-10AXL, SIL-20A, SIL-20AC, SIL-20ACHT, SIL-20ACXR, SIL-20AHT, SIL-20AXR, SIL-HTA, SIL-HTC	Shimadzu controller (see HPLC Other)	Ready
	Snir ^{/1,2}	Sion 120	RS232	Testing
	Spark ^{/4}	Alias	RS232	Ready
	Spark	Alias, Alias UHPLC (including User Program)	RS232, USB or LAN	Ready
	Spark	Basic Marathon (Plus, XT), Endurance, Midas, Triathlon	RS232	Ready
	Spark	Integrity	RS232, USB or LAN	Ready
	Sykam	S5200	RS232	Ready
	Sykam ^{/1}	S5250, S5300	RS232	Ready
	VICI - Valco	Multi-position valve used as AS	RS232 - UNI Ruby	Ready
New	Watrex ^{/4}	Streamline AS1, Streamline AS2	RS232	Ready
	YoungIn Chromass	ChroZen UHPLC ALS, YL2000H, YL2000HT, YL2100H, YL3000A, YL3100A, YL3200A	LAN	Ready
	YoungIn Chromass ^{/4}	YL9150	RS232	Ready
	YoungIn Chromass ^{/1}	YL9150 Plus	LAN	Testing
	YoungIn Chromass	YL9151	RS232	Ready
		Any autosampler with the Ready input and Inject output (Inject marker)	I/O device	Ready

HPLC Autosampler Accessories (p/n A26)

Producer	Name	Interface	Status
R-Biopharm	RIDACREST	RS232, USB or LAN	Ready
Spark	ACE, HPD, Symbiosis	RS232, USB or LAN	Ready

HPLC Column Ovens (p/n A24)

Producer	Name	Interface	Status	
Agilent	1100 TCC	GPIB, LAN	Ready	
Agilent	1200 TCC	LAN	Ready	
Agilent	G1316A, G1316B, G1316C	LAN - ICF	Ready	
Agilent	G7116A, G7116B, G7130A	LAN - ICF	Testing	
Analytical Technologies Limited ^{/4}	CO-3400U	RS232	Testing	
New	ANALYTIK Jena ^{/4}	PQ LC 4115	RS232	Testing
New	Chromsystems ^{/4}	CLC 360	RS232	Testing
	Cluzeau Info-Labo	Croco-CIL, Gecko-2000, Igloo-CIL	RS232 - UNI Ruby	Ready
New	CQS ^{/4}	Climax S 4115, Climax S 4120	RS232	Testing
	Ecom ^{/1}	ECO2080, ECO2099	RS232 or LAN	Ready
	Ecom	LCO 102	RS232 - UNI Ruby	Ready
New	Esensing	Mistral	RS232	Ready
New	Futecs ^{/1}	CT-6000	RS232 - UNI Ruby	Testing
	Hanon Instruments ^{/4}	LC7030	LAN	Ready
	Hitachi ^{/8,9}	Chromaster 5310, Chromaster Ultra RS 6310	USB	Ready
	Hitachi ^{/9}	L-2300	USB	Ready
	Hitachi ^{/6,8,9}	Primaide 1310	USB	Ready
	Knauer ^{/1}	CT2.1, PLATINblue T-1, S4000, S4050	RS232 or LAN	Ready

	PerkinElmer ^{/1}	Flexar LC Column Oven	RS232	Testing
	PerkinElmer	NexSAR Column Oven	RS232	Ready
	R-Biopharm	RIDACREST Column Oven	RS232	Ready
	Recipe	CH3000	RS232 - UNI Ruby	Ready
	Recipe	HT4000	RS232 or LAN	Ready
	RIGOL ^{/1}	L-3400	RS232	Testing
New	Schambeck ^{/4}	S 5120	RS232	Testing
	Shimadzu	CTO-10A, CTO-10AC, CTO-10ACvp, CTO-10ASvp, CTO-10Avp, CTO-20A, CTO-20AC	Shimadzu controller (see HPLC Other)	Ready
	Spark	Mistral	RS232	Ready
	SSI	605 COC	RS232 - UNI Ruby	Ready
New	Sykam ^{/1}	S4115, S4120	RS232	Testing
	Techlab	K-7	RS232 - UNI Ruby	Ready
	Thermotechnic Products ^{/4}	JetStream	RS232	Ready
	Unimicro ^{/1}	Unimicro Thermostat	RS232	Testing
	VICI - Valco ^{/1}	Fast Temperature Programmer	RS232	Ready
	Watrex	CTC 100	RS232 - UNI Ruby	Ready
New	Watrex ^{/1}	Streamline CT1	USB	Testing
	YoungIn Chromass ^{/1}	ChroZen HPLC-A CC, ChroZen UHPLC CC, YL9130	LAN	Ready

Fraction collectors (p/n A24)

	Producer	Name	Interface	Status
	Advantec	CHF122SC	RS232 - UNI Ruby	Ready
	Agilent	G1364A, G1364B, G1364C, G1364D, G1364E, G1364F, G5664A, G5664B, G7159B, G7166A	LAN - ICF	Testing
	Bio-Rad	BioFrac	A/D converter + UNI Ruby	Ready
	Büchi ^{/4,12}	C-660	RS232	Ready
	Ecom ^{/1}	ECF2096	RS232 or LAN	Testing
	Gilson	203B, 204	RS232 - UNI Ruby	Ready
Update	Japan Analytical Industry Co., LTD. ^{/1}	FC-7000	RS232	Ready
	Knauer ^{/1,12}	S3000, S3050	RS232 or LAN	Ready
	Labomatic Instruments ^{/4,12}	Labocol VARIO-4000, Labocol VARIO-4000 Plus	RS232 or LAN	Testing
	LKB Amersham	Frac100	A/D converter + UNI Ruby	Testing
	Rheodyne	MX Series II valves used as FC	USB	Ready
	Teledyne ISCO	FOXY 200, FOXY Jr., FOXY R1, FOXY R2	RS232 - UNI Ruby	Ready
	Teledyne ISCO ^{/4,12}	FOXY Jr., FOXY R1, FOXY R2		Ready
	VICI - Valco	Multi-position valve used as FC	RS232 - UNI Ruby	Ready
	Waters	Fraction Collector II	RS232 - UNI Ruby	Ready
	YoungIn Chromass	YL FC1, YL FC2	RS232 - UNI Ruby	Ready
		Any fraction collector using the Next and Collect/Waste events	Clarity-controlled device with 2 TTL outputs (INT7, etc.) + UNI Ruby	Ready

Valves (p/n A24 or A23)

	Producer	Name	Interface	Status
New	Agilent	G1156A, G1157A, G1158A, G1158B, G1159A, G1160A, G1162A, G1163A, G1170A, G4243A, G7170B, G9322A	LAN - ICF	Testing
New	ANALYTIK Jena ^{/4}	PQ LC 6510	RS232	Ready
New	Chromsystems ^{/4}	CLC 230	RS232	Ready
	Knauer ^{/1}	K-2, K-6, K-12, K-16, V2, V6, V12, V16, V2.1S, VU4.1	RS232 or LAN	Ready
	PerkinElmer	NexSAR Switching valve	RS232	Ready
	Rheodyne	LabPro Series valves	RS232 - UNI Ruby	Ready
	Rheodyne	MX Series II valves	USB	Ready
	Rheodyne	Titan Series valves	RS232	Ready
	Shimadzu	FCV-10AL, FCV-10ALvp, FCV-11AL, FCV-11ALS, FCV-15AL, LPGE Unit (all through selected LC-10/20 pumps only)	Shimadzu controller (see HPLC Other)	Ready
	Shimadzu	FCV-12AH, FCV-14AH (all through selected LC-20 column ovens only)	Shimadzu controller (see HPLC Other)	Ready
New	Sykam ^{/1}	S6510	RS232	Ready
	Upchurch	V-2500	RS232	Ready
	VICI - Valco	Any valve mounted on two-position microelectric actuator, multiposition microelectric actuator or universal actuator.	RS232	Ready

Amino Acid Analyzers

	Producer	Name	Interface	Status
	Ingos ^{/1}	AAA 500	LAN	Testing
	Sykam ^{/2}	S4300 Amino Acid Reaction Module System	RS232	Ready

Capillary Electrophoresis (p/n A31)

	Producer	Name	Interface	Status
	Prince ^{/1}	Next CE	LAN	Testing

Balances for Clarity EA Extension (p/n A30)

	Producer	Name	Interface	Status
	Mettler Toledo	AT, MT, UMT, Excellence	RS232	Ready
	Sartorius	CP, GC and GP series, Cubis series, M2P	RS232	Ready

Other (p/n A24 or A23)

Producer	Name	Interface	Status
DataApex	Column Usage Monitor	UNI Ruby	Ready

Interface types:

RS232	=	Serial communication. Specify the type of the cable or instrument.
LAN	=	Network communication. Specify the type of the cable or instrument.
USB	=	USB communication. Specify the type of the cable or instrument.
GPIB	=	hardware: Interface card for older Agilent 1100 system, Part No.: IAG11 (PCI interface) or IAG12 (USB interface)
GSIOC	=	hardware: Converter for Gilson serial communication, Part No.: IGLN1
Zebrick	=	hardware: External LC Control device (4x D/A & 4x D/F converter), Part No.: B30
OPT-USB	=	hardware: Optical to USB converter, Part No.: OPT01
A/D converter	=	hardware: A/D converter is needed to gather data because the instrument does not support digital acquisition. Part No.: U31, U32, U34
I/O device	=	hardware: Device capable of sending and receiving TTL signals governed by Clarity, for example any DataApex A/D converter. Part No.: U31, U32, U34
LC UNI	=	software: Universal LC control module for online pumps, Part No.: A24
FC GP	=	software: General Purpose Fraction Collector control module for Clarity, Part No.: A24
UNI RUBY	=	software: Universal control module for various instruments, Part No.: A24
ICF	=	software: Agilent ICF control Libraries embedded in Clarity, Part No.: A24, A26, A29 depending on the type of controlled instruments

Status:

Ready	=	the module has been subject to extensive testing and has been officially RELEASED.
Testing	=	the module is already functional and its development is in stage of external testing often called public BETA. It is possible to get this control module and help us with fine tuning the driver for your specific configuration.

Notes:

1. *Developed by instrument's producer.*
2. *Exclusively available from the instrument's producer in the particular OEM version.*
3. *Exclusively available from the instrument's producer – needs a release code.*
4. *Developed by third party.*
5. *Only models with Firmware 5.11 and lower.*
6. *Only limited number of these systems on one Clarity station is possible - see appropriate manual or product page.*
7. *Available from Dalian Elite company only (in the Chromsoft OEM version).*
8. *Exclusively available in selected territories through companies approved by Hitachi company.*
9. *One module in the system must have Interface Card installed.*
10. *The G1315A, G1315B, G1365A, G1365B sub-types only.*
11. *Selected models only - Check with Elysia-raytest company.*
12. *Available from Knauer company only (requires Knauer FRC option).*
13. *Only models with serial option - Check with Sedere company.*
14. *Only the detectors sold by Knauer company.*
15. *Exclusively available from the instrument's producer – only sold by instrument's producer.*